

LDI medical™

ON A MISSION TO REDUCE C. DIFF? Go with the #1 choice.

Introducing: EnviroSorb® CommodeShield™

You told the LDI Innovation Lab about the challenges of infection protection with patients using and staff cleaning commode buckets. Our EnviroSorb® CommodeShield™ has you covered as the only padded liner that completely covers the seat. This 3-ply padded liner removes moisture, eliminates splash back and prevents leaks.


- › INCREASES EFFICIENCY WITH EASY WRAP, CINCH & TOSS DESIGN
- › MAXIMIZE PATIENT COMFORT WITH SOFT SIT PADDING
- › LINER CAN HANDLE UP TO 1.1 LITER OF LIQUID
- › TURNS FLUIDS INTO GEL WITHIN SECONDS

Specifications

PART	ESCS-3524	ESCS-2824
SIZE	35" x 24"	28" x 24"
PACKAGING	5 each / bag, 10 bags / case	5 each / bag, 10 bags / case



BOOSTS INFECTION PROTECTION WITH 3 LAYERS:



- 1. Wicking Layer**
Removes moisture
- 2. Absorbent Layer**
Prevents splash back
- 3. Moisture Barrier Outer Layer:**
Prevents leaks