

STERILE

LDImedical™
WE LIVE FOR DEATHCARE™

WANT TO GET A HANDLE ON SAFETY? Get a glove with real grip.

Introducing: Surgeon Gloves - Chloroprene

You told the LDI Innovation Lab about the unique needs when it comes to protecting the hands of surgical professionals. Designed to be worn by staff in a sterile field, our Surgeon Gloves keep fluids away from the caregiver and protect against the transfer of microorganisms, particulates, etc from or to the caregiver in a clinical space.

- › POLYMER COATED INTERIOR FOR DAMP HAND DONNABILITY
- › IMPROVED GRIP WITH TEXTURED FINGERS

Specifications

PART	PPE-134(Size)
SIZE	5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
PACKAGING	25 pair / box, 4 boxes / case



Manufactured in compliance with:

- ASTM 6977
(Standard Specification for Polychloroprene Surgical Gloves)
- ASTM D3577
(Standard Specification for Rubber Surgical Gloves)
- ASTM D5151
(Standard Specification for Detection of Holes in Medical Gloves)
- ASTM F1671
(Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens)

Product Details: Polychloroprene material • 5.5 mil thick • White • Latex-free • Powder-free • Sterile

