


SURGICAL POSITIONERS NOT HOLDING UP? We're in a position to help.

Introducing: EnviroSoft® Surgical Positioners

You told the LDI Innovation Lab about the challenges with existing surgical positioning devices. That's why we created Envirosoft® Surgical Positioners - the soft, durable, and reusable option that offers maximum support for the patient and cost savings for the facility ... without the mess (or duct tape).

- › MAXIMIZE SUPPORT OF SOFT TISSUE WITH MINIMAL PRESSURE ON SKIN
- › DOESN'T LEAK AND USEABLE EVEN WHEN NICKED OR CUT
- › SPEND LESS MONEY WITH A COST-EFFECTIVE, REUSABLE OPTION
- › EASY TO CLEAN SIMPLY WIPE DOWN BETWEEN PROCEDURES



Specifications

	CHEST ROLL	HEEL CUP	SUPINE HEAD BLOCK	ARM CRADLE	PRONE HEAD BLOCK
PART	ESP-CP	ESP-HC	ESP-HPS	ESP-AP	ESP-PP
SIZE	16"L x 4"W x 3"H	5 3/4"L x 1 3/4"H	7 7/8"L x 3"H	16"L x 5.5"W x 2" H	9.5"W x 6"H
PACKAGING	10 each / case	6 each / case	30 each / case	10 each / case	5 each / case

Product Details: Blue • TPE construction • PVC-free • Latex -free

