

FLUIDS ON THE FLOOR ARE JUST BAD... So we created this pad.

Introducing: EnviroSorb® Original Surgical Floor Pads

You told the LDI Innovation Lab that fluid-intensive surgical procedures are hard enough without the mess. Our EnviroSorb® Original Surgical Floor Pads let you focus on tissue recovery or surgery by providing a super-absorbent barrier that keeps fluids contained—and you from slipping.

- › ABSORBS UP TO 1 GALLON PER SQ. YARD
- › KEEPS THE FLOOR CLEAN AND DRY WITH PROTECTIVE BACKING
- › HELPS YOU AVOID SLIPS AND FALLS DURING TISSUE RECOVERY

Specifications

PART	ESP-1616	ESP-3630	ESP-3640	ESP-3072	ESP-4672
SIZE	16" x 16"	36" x 30"	36" x 40"	30" x 72"	46" x 72"
PACKAGING	50 pads / case	25 pads / case	25 pads / case	25 pads / case	10 pads / case

