

STRUGGLING WITH URINE ODORS.

You're in luck.

We've got a pad for that.

Introducing: EnviroSorb® Washable Underpads with Urine Odor Control

You told the LDI Innovation Lab about the challenges of controlling odor and maintaining patient dignity. We responded by developing EnviroSorb® Washable Underpads with Urine Odor Control - the simple and cost-effective solution that reduces the smell of urine incontinence, while keeping patients comfortable and dignified.



- › MINIMIZES URINE ODORS WITH PROPRIETARY URINE-NEUTRALIZING MATERIAL
- › REDUCES RISK OF BEDSORES & DECUBITUS ULCERS BY WICKING FLUIDS AWAY FROM PATIENT SKIN
- › CAN BE LAUNDERED WITH EITHER BLEACH OR PEROXIDE

Specifications

PART	UOC-3436
SIZE	34" x 36"
PACKAGING	24 each / case

Product Details:

- Quilted and soft, but heavyweight, 55/45 cotton/polyester twill face with wicking properties
- Polyester/rayon soaker core treated with urine odor neutralizer
- Green fluid barrier backing

