

LOOKING TO PROTECT PATIENTS AND STAFF? Let us lend a hand.

Introducing: Exam Gloves

You told the LDI Innovation Lab about the risks inherent with fluids during procedures. Worn by staff to cover their hands, our Exam Gloves keep fluids away from the caregiver and protect against the transfer of microorganisms, particulates, etc from or to the caregiver in a clinical space.

- › IMPROVED GRIP WITH TEXTURED FINGERS
- › EXCEPTIONAL COMBINATION OF COMFORT, PROTECTION, AND VALUE
- › ENHANCED DORNING STRENGTH WITH BEADED CUFF



Specifications

PART	NEG-PF3-(Size)
SIZE	S, M, L, XL
PACKAGING	100 each / box, 10 boxes / case

Product Details: Nitrile material • 3 mil thick • Violet Blue • Latex-free • Powder-free • Non-sterile

