


A PAD THAT NEUTRALIZES FORMALDEHYDE DURING TISSUE DISSECTION? Let that sink in for a bit.

Introducing: EnviroSorb® Formaldehyde Neutralizing Pads

You told the LDI Innovation Lab about the risks associated with toxic formaldehyde. And dissection happens. That's why we created EnviroSorb® Formaldehyde Neutralizing Pads. Made with a treated soaker fabric that chemically neutralizes formaldehyde, these pads reduce your exposure during dissection and make cleanup far easier.



- › ABSORBS & CONTAINS FORMALDEHYDE FLUID FROM SPILLS & DURING DISSECTION
- › REDUCES EXPOSURE TO HARMFUL VAPOR BY CHEMICALLY NEUTRALIZING FORMALDEHYDE
- › IMPROVES WORKPLACE SAFETY & OSHA COMPLIANCE

Specifications

PART	FNP-0810B-400
SIZE	8" x 10"
PACKAGING	400 each / case

Product Details: White soaker fabric • Blue scrim top layer

