


# A PAD THAT NEUTRALIZES FORMALDEHYDE DURING TISSUE DISSECTION? Let that sink in for a bit.

## Introducing: EnviroSorb® Formaldehyde Neutralizing Pads

You told the LDI Innovation Lab about the risks associated with toxic formaldehyde. And dissection happens. That's why we created EnviroSorb® Formaldehyde Neutralizing Pads. Made with a treated soaker fabric that chemically neutralizes formaldehyde, these pads reduce your exposure during dissection and make cleanup far easier.



- › ABSORBS & CONTAINS FORMALDEHYDE FLUID FROM SPILLS & DURING DISSECTION
- › REDUCES EXPOSURE TO HARMFUL VAPOR BY CHEMICALLY NEUTRALIZING FORMALDEHYDE
- › IMPROVES WORKPLACE SAFETY & OSHA COMPLIANCE

## Specifications

<b>PART</b>	FNP-0810B-400	FNP-1620B-100
<b>SIZE</b>	8" x 10"	16" x 20"
<b>PACKAGING</b>	400 each / case	100 each / case

**Product Details:** White soaker fabric • Blue scrim top layer

