

# Autopsy Twine/ Ligature.

Breakage-free.  
Cleaner. Safer. Smarter.



## A STRONGER SOLUTION FOR POST-MORTEM SUTURING

Introducing the Autopsy Twine/Ligature. A durable, waxed solution for tissue recovery, autopsy, and other post-mortem procedures, specifically designed to:

- Easily thread through the eye of the needle.
- Resist breaking during post-mortem suturing.
- Maintain position due to its waxed and fluid-resistant construction.
- Ensure safety with a smoother sew.



Built to help you maintain position during suturing, the Autopsy Twine/Ligature keeps your team safe and keeps every procedure clean and dignified.

## PRODUCT DETAILS

5 CORD LINEN MATERIAL  
WAXED FINISH  
NATURAL COLOR  
APPROX. 850 YARDS

## SPECIFICATIONS

PART	AT-5W-LINEN
SIZE	1 LB. SPOOL
PACKAGING	PACKAGED INDIVIDUALLY