

A BODY BAG THAT ABSORBS FLUIDS? Soak that up for a minute.

INTRODUCING:

EnviroSorb® Body Bag

You told the LDI Innovation Lab about the nightmare of opening body bags full of fluids. That's why we created the EnviroSorb® Body Bag. This durable, PVC-Free, bag has a built-in liner that can absorb 1.5 gallons of fluid. Building on the success of our EnviroMed-Bag®, it provides even more exposure protection while eliminating cleanup during removal, transport, and storage of a body.

- › BUILT-IN LINER ABSORBS 1.5 GALLONS OF FLUID
- › PREVENTS LEAKS AND EXPOSURE TO BLOOD-BORN PATHOGENS
- › SAFE FOR CREMATION

Specifications

PART	EMD2-ESP-3692A
SIZE	36" x 92"
PACKAGING	6 bags / case

Product Details: White • 8 mil thick • 3 tags included • Chlorine Free • PVC-Free

