

# EnviroSorb Body Bag.

Built-In Fluid Absorption.  
Cleaner. Safer. Smarter.



## A SMARTER SOLUTION FOR FLUID MANAGEMENT DURING TRANSPORTATION

Introducing the EnviroSorb Body Bag. A durable, PVC-free and chlorine-free solution designed to:

- Absorb fluids with a built-in liner holding up to 1.5 gallons.
- Prevent leaks and exposure to blood-borne pathogens.
- Simplify transport and storage by eliminating messy cleanup.
- Ensure safety with materials safe for cremation.



Built on the success of our EnviroMed Body Bag and EnviroSorb Pads, Liners, and Wraps, the EnviroSorb Body Bag keeps your team safe and every transport clean and dignified.

### PRODUCT DETAILS

WHITE  
8 MIL THICK  
3 TAGS INCLUDED  
CHLORINE & PVC FREE

### SPECIFICATIONS

PART	EMD2-ESP-3692A
SIZE	36" X 92"
PACKAGING	6 BAGS/CASE