

EnviroSorb CommodeShield™

Infection Protection.
Cleaner. Safer. Smarter.



A SMARTER SOLUTION FOR WASTE MANAGEMENT

Introducing the EnviroSorb CommodeShield™. Our unique design boosts infection protection and is designed to:

- Ensure safety with a 3-ply padded liner.
- Hold up to 1.1 liters of liquid; turn fluids into gel within seconds.
- Maximize patient comfort with soft-sit padding.
- Simplify cleanup with easy wrap, cinch, and toss design.



Our EnviroSorb CommodeShield™ offers 3 layers of protection (wicking layer removes moisture, absorbant layer prevents splash, and moisture barrier layer prevents leaks) to keep your team safe and maintain patient dignity.

PRODUCT DETAILS

3-PLY PADDED LINER
COMPLETELY COVERS SEAT

SPECIFICATIONS

PART	ESCS-3524	ESCS-2824
SIZE	35" X 24"	28" X 24"
PACKAGING	PACK OF 5, CASE OF 50	PACK OF 5, CASE OF 50