


EnviroSorb® Original Floor Pad.

Mess-Free Fluid Absorption.
Cleaner. Safer. Smarter.



AN EFFECTIVE SOLUTION FOR FLUID MANAGEMENT DURING RECOVERY

Introducing the EnviroSorb® Floor Pad. A durable pad designed to:

- Absorb up to 1 gallon of fluid per square yard.
- Keep the floor clean and dry with a protective backing.
- Eliminate messy cleanup.
- Ensure safety by preventing slips and falls during tissue recovery.



Developed to keep your team safe and keep every recovery clean and dignified.

PRODUCT DETAILS

WHITE
NON-SKID BACKING

SPECIFICATIONS

PART	ESP-3630	ESP-3640	ESP-3072	ESP-4672
SIZE	36" X 30"	36" X 40"	30" X 72"	46" X 72"
PACKAGING	25 PADS/CASE	25 PADS/CASE	25 PADS/CASE	10 PADS/CASE