


EnviroSorb Super Capacity Floor Pad.

Volume Fluid Absorption.
Cleaner. Safer. Smarter.



AN EFFECTIVE SOLUTION FOR FLUID MANAGEMENT DURING RECOVERY

Introducing the EnviroSorb Super Capacity Floor Pad. A heavy-duty pad designed to:

- Absorb and solidify up to 2 gallons of fluid per square yard.
- Keep the floor clean and dry with a protective backing.
- Eliminate messy cleanup.
- Ensure safety by preventing slips and falls during tissue recovery.



Built on the success of our EnviroSorb Original Floor Pad, the Super Capacity Floor Pad keeps your team safe and keeps every recovery clean and dignified.

PRODUCT DETAILS

WHITE
NON-SKID BACKING

SPECIFICATIONS

PART	ESP-SAP-3640	ESP-SAP-3072
SIZE	36" X 40"	30" X 72"
PACKAGING	25 PADS/CASE	25 PADS/CASE