

# Post-Mortem Head Block.

Maintaining Dignity.  
Cleaner. Safer. Smarter.



## A SMARTER SOLUTION FOR FLUID MANAGEMENT DURING PROCEDURES

Introducing the Post-Mortem Head Block. A single-use solution designed to:

- Raise the head to reduce pooling of blood and fluids in the eyes and facial capillaries.
- Create a more stable setting for work.



The Post-Mortem Head Block ensures that every decedent is treated with care and respect.

### PRODUCT DETAILS

WHITE  
POLYSTYRENE MATERIAL  
SINGLE USE

### SPECIFICATIONS

PART	DHB-01
SIZE	UNIVERSAL
PACKAGING	24 EACH/CASE