

STERILE

LDImedical™
WE LIVE FOR DEATHCARE™

SURGERY IS HARD ON THE HANDS. Get a softer glove.


Introducing: Surgeon Gloves - Polyisoprene

You told the LDI Innovation Lab about the unique needs of protecting the hands of surgical professionals from the transfer of microorganisms in a clinical space. Designed to be worn by staff in a sterile field these gloves are made from extra durable polyisoprene and most closely mimic the soft, supple, comfortable fit and feel of a latex glove, without the problems associated with latex sensitivity.

- › UNIQUE HYDROGEL INTERIOR FOR DAMP HAND DONNABILITY
- › IMPROVED GRIP WITH TEXTURED FINGERS

Specifications

PART	PPE-146(Size)
SIZE	5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
PACKAGING	50 pair / box, 4 boxes / case



Manufactured in compliance with:

- ASTM D3577
(Standard Specification for Rubber Surgical Gloves)
- ASTM D5151
(Standard Specification for Detection of Holes in Medical Gloves)
- ASTM F1671
(Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens)

Product Details: Polyisoprene material • 9.8 mil thick • White • Latex-free • Powder-free • Sterile

