

STERILE

LDImedical™
WE LIVE FOR DEATHCARE™

NEED EXTRA FLUID PROTECTION? We've got you covered.

Introducing: Surgical Gowns - Level 3 Reinforced

You told the LDI Innovation Lab about the enhanced protection required by your surgical staff. Our reinforced Surgical Gowns are worn by staff to cover their arms and torso, down to below the knees during surgical procedures in a sterile field. These Level 3 Reinforced gowns keep fluids away from the caregiver and protect against the transfer of microorganisms, particulates, etc from or to the caregiver in a clinical space.

- > COMPLIES WITH AAMI PB 70 STANDARD (SURGICAL DRAPES & GOWNS)
- > MEETS AAMI LEVEL III FOR WATER RESISTANCE
- > BALANCE COMFORT AND PROTECTION WITH POLYETHYLENE REINFORCEMENT IN HIGH SPLASH AREAS
- > TOWEL INCLUDED



Product Details:

- SMS polypropylene material
- Polyethylene coating
- Blue
- Knit cuffs
- Latex-free
- Sterile
- Meets AAMI Level III for water resistance

Specifications

PART	PPE-E100350	PPE-E100360	PPE-E100370
SIZE	L	XL	XXL
PACKAGING	36 each / case	32 each / case	28 each / case

