

STERILE

LDImedical™
WE LIVE FOR DEATHCARE™

NEED MORE PROTECTION? Get the gown that gets it done.

Introducing: Surgical Gowns - Level 4

You told the LDI Innovation Lab about the enhanced protection required by your surgical staff. Our reinforced Surgical Gowns are worn by staff to cover their arms and torso, down to below the knees during surgical procedures in a sterile field. These Level 4 gowns keep fluids away from the caregiver and protect against the transfer of microorganisms, particulates, etc from or to the caregiver in a clinical space.

- › COMPLIES WITH AAMI PB 70 STANDARD (SURGICAL DRAPES & GOWNS)
- › MEETS AAMI LEVEL IV FOR SYNTHETIC BLOOD AND VIRAL PENETRATION
- › TOWEL INCLUDED



Product Details:

- SMS polypropylene material
- Polyethylene coating
- Blue
- Knit cuffs
- Latex-free
- Sterile
- Meets AAMI Level IV for synthetic blood and viral penetration

Specifications

PART	PPE-E104050	PPE-E104060	PPE-E104070
SIZE	L	XL	XXL
PACKAGING	44 each / case	44 each / case	32 each / case

