



FINALLY STOPPING THE LEAKS... It's in the bag.

Introducing: EnviroMed-Bag® Body Bag

You told the LDI Innovation Lab you were tired of ripping and leaking body bags. That's why we created the EnviroMed-Bag® lineup. Designed with super-durable PVC-Free materials and heat-sealed seams, our bags reduce tears and fluid leakage during removal, transport and storage of a body. The end result is a product that minimizes cleanup and protects from exposure to hazardous fluids.

- › AVOIDS THE MESS WITH LEAK & TEAR RESISTANT DESIGN
- › PREVENTS LEAKS AND EXPOSURE TO BLOOD-BORN PATHOGENS
- › SAFE FOR CREMATION



Specifications

	INFANT	PEDIATRIC	ADULT	XL ADULT*	BARIATRIC
PART	EMD2-1428A	EMD2-2448A	EMD2-3692A	EMD2-4096A	EMD2-8296A
SIZE	14" x 28"	24" x 48"	36" x 92"	40" x 96"	82" x 96"
PACKAGING	12 bags / case	12 bags / case	12 bags / case	12 bags / case	6 bags / case

Product Details: White • 8 mil thick • 3 tags included • Chlorine Free • PVC-Free

*XL ADULT also available w/straps + wraps: EMD2-4096M