

PROFESSIONAL APPEARANCE. BETTER GRIP. You have to hand it to us.

INTRODUCING:

MICROFLEX® Black Dragon® Zero Nitrile Exam Gloves

At LDI Medical, we recognize the importance of medical professionals presenting a calm, professional appearance during sensitive procedures, especially when working with decedents and their families. The MICROFLEX® Black Dragon® Zero is a black medical exam glove designed to mask stains during use. With textured fingertips, these gloves help provide a secure grip and are designed with a soft, durable nitrile formulation for elevated comfort and strength.

- › BLACK NITRILE MASKS STAINS & MAINTAINS PROFESSIONAL APPEARANCE
- › TEXTURED FINGERTIPS PROVIDE A SECURE GRIP
- › SOFT FORMULATION FOR GREAT COMFORT & DURABILITY



Specifications

PART	BD-1001-NPF	BD-1002-NPF	BD-1003-NPF	BD-1004-NPF	BD-1005-NPF
SIZE	Small	Medium	Large	XL	2XL
PACKAGING	100 gloves/box, 10 boxes/case	100 gloves/box, 10 boxes/case	100 gloves/box, 10 boxes/case	100 gloves/box, 10 boxes/case	100 gloves/box, 10 boxes/case

Product Details: Black • 3.5 mils palm thickness • 5.5 mils finger thickness