

LDI medical™

# COMFORT MEETS DURABILITY.

## These gloves get the job done.

**INTRODUCING:**

### MICROFLEX® Cobalt® Nitrile Exam Gloves

You told the LDI Innovation Lab about the challenges of finding the best day-to-day glove. These MICROFLEX® Cobalt® Exam Gloves offer consistent, dependable comfort and durability. Fully textured, these tough, 4 mil nitrile gloves help wearers securely grip tools and materials. These are medical exam graded and tested for use with chemotherapy drugs, offering enhanced protection that provides a good option for healthcare and deathcare professionals alike. These gloves also provide protection from Type I latex allergy.

- › 4 MIL COBALT NITRILE
- › FULLY TEXTURED FOR AN ENHANCED GRIP
- › PROTECTION FROM TYPE I LATEX ALLERGY



## Specifications

PART	N191	N192	N193	N194
SIZE	Small	Medium	Large	XL
PACKAGING	100 gloves/box, 10 box/case	100 gloves/box, 10 box/case	100 gloves/box, 10 box/case	100 gloves/box, 10 box/case

**Product Details:** Cobalt • 3.9 mils palm thickness • 4.3 mils finger thickness