

LDI medical™

EMBALMING & AUTOPSY ARE TOUGH. Get a glove that can handle it.

INTRODUCING:

MICROFLEX® SafeGrip™ Latex Exam Gloves

You told the LDI Innovation Lab about the challenges of damaged gloves. These MICROFLEX® SafeGrip™ latex medical exam gloves are 3x thicker than standard disposable latex gloves, ensuring greater resistance against snags and tears. Perfect for protection during heavy-duty applications, these thick latex gloves' soft, elastic formulation delivers comfort and flexibility. Gloves feature textured fingertips with enhanced grip performance. Also, donning is made easy thanks to a double-chlorinated design.

- › 3X THICKER, GREATER SNAG AND TEAR RESISTANCE
- › SOFT LATEX FORMULATION: A MORE ELASTIC FIT & SUPERIOR COMFORT
- › WITH AN EXTENDED CUFF, THESE GLOVES OFFER ENHANCED FOREARM & WRIST PROTECTION



Specifications

PART	SG-375-S	SG-375-M	SG-375-L	SG-375-XL
SIZE	Small	Medium	Large	XL
PACKAGING	50 gloves/box, 10 box/case	50 gloves/box, 10 box/case	50 gloves/box, 10 box/case	50 gloves/box, 10 box/case

Product Details: Blue • 11 mils palm thickness • 14.2 mils finger thickness